



## Product Specification

**Product Information:**

<b>Products:</b>	Sulfuric acid 0.02N
<b>Catalog:</b>	VS17252
<b>Grade:</b>	Standard volumetric solution

**Reagent Used:**

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Sulfuric acid	7664-93-9	ACS

**Specification:**

Test	Specification
Appearance	Colorless liquid
Assay at 20°C	0.02 N $\pm$ 5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

