



## Product Specification

### Product Information:

**Products:** Sulfuric acid 0.01N (0.005M)

**Catalog:** VS17264

**Grade:** Standard volumetric solution

### Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Sulfuric acid	7664-93-9	ACS

### Specification:

Test	Specification
Appearance	Clear colorless solution
Assay	0.01 N $\pm$ 5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

