



Product Specification

Product Information:

Products:	Zinc sulfate 1M (2N)
Catalog:	VS17277
Grade:	Standard volumetric solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Zinc sulfate heptahydrate	7446-20-0	ACS

Specification:

Test	Specification
Appearance	Clear colorless liquid
Assay at 20°C	1.0 M \pm 5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

