



Product Specification

Product Information:

Products:	Iodine solution 0.05 N
Catalog:	VS17287
Grade:	Standard volumetric solution

Reagent Used:

Reagents	CAS #	Grade
Water	7732-18-5	USP/ASTM
Iodine	7553-56-2	ACS
Potassium iodide	7681-11-0	ACS
Hydrochloric acid	7647-01-0	ACS

Specification:

Test	Specification
Appearance	Dark red brown solution
Assay	0.05 N \pm 5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

