



# Product Specification

Product : Water  
Catalog# : LC9883  
CAS : 7732-18-5  
Formula : H<sub>2</sub>O  
MW : 18.02  
Grade : HPLC

TEST	SPECIFICATION
Description	Colorless liquid.
Residue on evaporation	≤ 1 ppm
Specific conductance @ 25°C	≤ 1.5 μS/cm
<b>UV absorbance</b>	
@ 190 nm	≤ 0.01
@ 200 nm	≤ 0.01
@ 210 nm	≤ 0.01
@ 250 nm	≤ 0.01
@ 400 nm	≤ 0.01

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

