



Product Specification

Product : Bromocresol Green Free Acid
Catalog# : I10503
CAS : 76-60-8
Formula : $C_{21}H_{14}Br_4O_5S$
MW : 698.04

TEST	SPECIFICATION
Appearance	Beige yellow to yellow brown powder
Loss on drying	< 3%
Clarity of Solution	Passes Test
Aqueous solubility	To pass test
Visual Transition Interval	From pH 3.8 (Greenish yellow to yellow) To pH 5.4 (Blue)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

