



# Product Specification

Product : Bromocresol Green Sodium Salt  
Catalog# : I10504  
CAS : 62625-32-5  
Formula : C<sub>21</sub>H<sub>13</sub>Br<sub>4</sub>NaO<sub>5</sub>S  
MW : 720.02

TEST	SPECIFICATION
Description	Dark-brown to dark red or Greenish-black crystalline powder
Clarity of solution	To pass test
Absorptivity 1%, 1cm, in Water at $\lambda$ max. 616.2 nm	> 350
Visual Transition Interval	From pH 3.8 (yellow) To pH 5.4 (blue)
Loss on drying	< 5%
Solubility	H <sub>2</sub> O 40mg/ml
	EGME 50mg/ml
	ETOH 60mg//ml

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

