



Product Specification

Product : Bromophenol Red Sodium Salt
Catalog# : I10513
CAS : 102185-50-2
Formula : C₁₉H₁₁Br₂O₅SNa
MW : 534.15

TEST	SPECIFICATION
Appearance	Red to Redish-Brown powder
Clarity of Solution	Clear to a Slight Haze
Loss on drying (at 110°C /1 hr.)	< 5.0%
Visual Transition Interval	pH 5.2 (yellow-red) to pH 7.0 (purple)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

