



Tri Stains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Cresol Red Free acid
Catalog# : I10526
CAS : 1733-12-6
Formula : $C_{21}H_{18}O_5S$
MW : 382.43

TEST	SPECIFICATION
Appearance	Green to Maroon-Red to Dark Brown Powder
Assay	$\geq 95 \%$
Visual Transition Intervals	pH 0.2 (pink) to pH 2.0 (yellow) pH 7.2 (yellow) to pH 8.8 (red-purple)
Loss on Drying	$\leq 3 \%$
Solubility 0.1% (95% EtOH)	Clear yellow solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

