



Product Specification

Product : Methyl Red Sodium Salt
Catalog# : I10541
CAS : 845-10-3
Formula : $C_{15}H_{14}N_3O_2Na$
MW : 291.29

TEST	SPECIFICATION
Description	Dark Red to brownish yellow to brown powder
Melting Point	179-182° C
Clarity of alcohol Solution	To Pass test
Clarity of aqueous Solution	To Pass Test
Visual Transition Interval	From pH 4.2 (pink) To pH 6.2 (yellow)
Loss on drying @110°C /1hr	< 5%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

