



Product Specification

Product : Methylthymol Blue sodium salt
Catalog# : I10542
CAS : 1945-77-3
Formula : $C_{37}H_{40}N_2O_{13}Na_4S$
MW : 844.74

TEST	SPECIFICATION
Appearance	Brown to Dark-Green to black powder
Suitability as a metal indicator	To Pass Test
Clarity of solution	To Pass Test
Visual transition interval	From pH 6.5 (yellow) To pH 8.5 (blue)
Loss on drying @110°C	≤ 10 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

