



Product Specification

Product : Xylenol Blue
Catalog# : I10574
CAS : 125-31-5
Formula : C₂₃H₂₂O₅S
MW : 410.49

TEST	SPECIFICATION
Appearance	Dark-brown to purple powder
Melting Point	212 °C
Dye content	~90%
Visual Transition Interval(acid)	From pH 1.2 (red) To pH 2.8 (yellow)
Visual Transition Interval(alkaline)	To pH 9.8 (violet)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

