



Tri Stains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Xylenol Orange Tetrasodium Salt
Catalog# : I10575
CAS : 3618-43-7
Formula : $C_{31}H_{28}Na_5O_{13}SN_2$
MW : 760.58

TEST	SPECIFICATION
Appearance	Dark red to dark brown crystalline powder
Loss on drying (at 110 °C/1hr)	≤ 10.0 %
Clarity of solution (0.1 g/100 ml of water)	Clear Solution
Suitability for zinc Titration	To Pass Test
Absorptivity 1%/1cm, at λ max. 576.6 nm	>550

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

Tri Stains

