



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Amido Black 10B
Catalog# : S10023
CAS : 1604-48-8
CI number : 20470
Formula : C₂₂H₁₄N₆Na₂O₉S₂
MW : 616.50

TEST	SPECIFICATION
Appearance	Dark brown to reddish-black powder
Identification	To conform to structure and standard
Absorptivity 1%/1cm at Lambda Max 618.0 nm	750 - 960
Solubility 0.1% (dist.water)	Clear blue solution
Loss on drying @ 105 deg. C	10.0% maximum

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

