



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Azure B (Azure 1)
Catalog# : S10028
CAS : 531-55-5
CI number : 52010
Formula : C₁₅H₁₆ClN₃S
MW : 305.83

TEST	SPECIFICATION
Appearance	Brown to Dark green or Black powder
Assay	≥ 75 %
Solubility (0.1% Aq)	Clear blue solution
Absortivity 1% 1cm At Lambda max. 645nm	> 1500
Loss on drying at 110 °C (1Hr.)	≤ 8 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

