



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Azur II Eosinate
Catalog# : S10031
CAS : 53092-85-6
CI number : -
Formula : -
MW : -

TEST	SPECIFICATION
Product comments	A mixture of Azur I, Methylene Blue and Eosin Y
Appearance	Dark green to green-black powder
Infrared Spectrum	Conforms to structure and standard
Lambda max.	647 (524) nm in methanol
Solubility: (0.1% Acetic Acid)	Opaque, dark blue solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

