



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Biebrich Scarlet
Catalog# : S10033
CAS : 4196-99-0
CI number : 26905
Formula : C₂₂H₁₄N₄O₇S₂Na₂
MW : 556.49

TEST	SPECIFICATION
Description	Deep-red to dark brown powder
Identification	To conform to standard
Lambda Max.	505 nm in Methanol
Assay	≥ 90 %
Loss on drying: (110° C for 1 Hr.)	< 15%
Solubility (0.1% dist. H ₂ O)	Clear red solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

