



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Bismark Brown R
Catalog# : S10034
CAS : 5421-66-9
CI number : 21010
Formula : C₂₁H₂₄N₈
MW : 461.40

TEST	SPECIFICATION
Appearance	Brownish-black powder
Lambda max.	468nm in 50% ethanol
Solubility	To pass test
Loss on drying at 60 °C	≤ 1.0% w/w
Melting Point	222 °C

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

