



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Brilliant Black BN  
Catalog# : S10036  
CAS : 2519-30-4  
CI number : 28440  
Formula : C<sub>28</sub>H<sub>17</sub>N<sub>5</sub>Na<sub>4</sub>O<sub>14</sub>S<sub>4</sub>  
MW : 867.69

TEST	SPECIFICATION
Appearance	Black powder
Dye Content	≥ 60%
Water insoluble matter	≤ 0.2%
Un sulfonated Aromatic Amines	≤ 0.01%
Lead (Pb)	≤ 2.00 mg/kg
Mercury (Hg)	≤ 1.00 mg/kg
Heavy metals	≤ 40.00 mg/kg
Organic compound other than Coloring matter, percent by mass, Max.	≤ 0.80%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

