



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Brilliant Blue R250
Catalog# : S10038
CAS : 6104-59-2
CI number : 42660
Formula : C₄₅H₄₄N₃NaO₇S₂
MW : 825.99

TEST	SPECIFICATION
Appearance	Dark-Purplish-Brown Powder
Loss on Drying	≤ 5%
Dye content	≥ 80%
Solubility (10 mg/ml of water)	Dark-blue-black solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

