



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Brilliant Yellow
Catalog# : S10040
CAS : 3051-11-4
CI number : 24890
Formula : C₂₆H₁₈N₄Na₂O₈S₂
MW : 624.55

TEST	SPECIFICATION
Appearance	Brick red to rust colored powder
Dye content	~85%
Visual Transition Interval	From pH 6.4 (yellow) To pH 8.0 (red-orange)
Solubility	Soluble in water (2.5 gm/100 ml); slightly soluble in acetone and cellosolve.

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

