



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Chromotrope 2R
Catalog# : S10051
CAS : 4197-07-3
CI number : 16570
Formula : $C_{16}H_{10}N_2Na_2O_8S_2$
MW : 468.37

TEST	SPECIFICATION
Description	Brownish-red powder
Assay	>75%
Absorptivity 1%/1cm, (in Methanol) at λ max. 511.6 nm	>485
Loss on drying	< 10%
Solubility: 0.1% in water	Clear red solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

