



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Chromotrope FB  
Catalog# : S10052  
CAS : 3567-69-9  
CI number : 14720  
Formula :  $C_{20}H_{12}N_2Na_2O_7S_2$   
MW : 502.43

TEST	SPECIFICATION
Appearance	Dark Red to Brown to Black powder
Infrared Spectrum	Conforms to Structure
UV-Visible Spectrum	0.03 G/L, H <sub>2</sub> O
	E (516 +/- 3NM) = >14,500
	E (323 +/- 3NM) = >9,000
	E (217 +/- 3NM) = >33,000

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

