



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Cresyl Violet Acetate Certified
Catalog# : S10056
CAS : 10510-54-0
CI number : -
Formula : $C_{18}H_{15}N_3O_3$
MW : 321.34

TEST	SPECIFICATION
Description	Dark-green to Dark Green-black to Black powder
Identification	To pass test
Lambda Max	596 - 599 nm (in 50 % Ethanol)
Dye content	> 65.0%
Melting Point	140 – 145° C
Solubility	2% in water 1% in ethanol

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

