



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Crocein Orange G
Catalog# : S10057
CAS : 1934-20-9
CI number : 15970
Formula : $C_{16}H_{11}N_2NaO_4S$
MW : 350.32

TEST	SPECIFICATION
Description	Red-orange to orange powder
Identification	To pass test
Loss on drying @110 deg. C/1hr	< 20.0%
Solubility	Soluble in water; moderately soluble in acetone, methyl Cellosolve and ethanol
Absorptivity 1%/1cm (Water) At Lambda Max	485.0 nm in water

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

