



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Crystal Scarlet
Catalog# : S10059
CAS : 2766-77-0
CI number : 16250
Formula : $C_{20}H_{12}N_2Na_2O_7S_2$
MW : 502.43

TEST	SPECIFICATION
Description	Dark brown powder
Identification	To pass test
Melting Point	> 300 ° C
Lambda Max	509-514 nm
Dye content	~ 80%
Solubility 0.1% (Water)	Clear red solution
Loss on drying: at 110 Deg. C/ 1 hour	<20%

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

