



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Eriochrome Blue Black R
Catalog# : S10072
CAS : 2538-85-4
CI number : 15705
Formula : $C_{20}H_{13}N_2O_5SNa$
MW : 416.39

TEST	SPECIFICATION
Identification	To conform to standard
Solubility in Water / Ethanol	Passes Test
Specific Extinction (C-20mg/L)	
E @ 215nm	26,600 min.
E @ 513nm	7,660 min.
Dye Content	≥ 70 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

