



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Ethyl Eosin Certifiable
Catalog# : S10077
CAS : 6359-05-3
CI number : 45386
Formula : $C_{22}H_{11}Br_4KO_5$
MW : 714.07

TEST	SPECIFICATION
Appearance	Dark-red powder
Lambda max.	532nm in 50% ethanol
Loss on drying: (at 105 Deg. C)	≤ 5 %
Solubility: (1 mg/ml 95% ETOH)	Clear pink-orange solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

