



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Ethyl Violet
Catalog# : S10078
CAS : 2390-59-2
CI number : 42600
Formula : $C_{31}H_{42}N_3Cl$
MW : 492.15

TEST	SPECIFICATION
Appearance	Olive-Green Powder
Visual Transition Interval	pH 0.0 (yellow) to pH 3.6 (blue)
Solubility	Slightly soluble in water Soluble in methanol and ethanol
Dye content	≥ 80 %
Insolubles	≤ 1.0 %
Moisture	≤ 4.0 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

