



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Fast Garnet GBC Base
Catalog# : S10080
CAS : 97-56-3
CI number : 11160
Formula : $C_{14}H_{15}N_3$
MW : 225.29

TEST	SPECIFICATION
Appearance	Orange to red-brown to dark brown powder
Assay	>97%
Identification (IR)	To conform to structure
Loss on drying	<3%
Melting Point	100-103° C

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

