



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Fast Sulphon Black F  
Catalog# : S10082  
CAS : 3682-47-1  
CI number : 26990  
Formula :  $C_{30}H_{17}N_4O_{11}S_3Na_3$   
MW : 774.6

TEST	SPECIFICATION
Appearance	Black powder
Purity (TLC)	≥ 95.0 %
Application	Suitable as an indicator for EDTA titration
Solubility	Dark green solution at 4 mg plus 4.0 ml Of water

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

