



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Fat Brown RR
Catalog# : S10083
CAS : 6416-57-5
CI number : 11285
Formula : $C_{16}H_{14}N_4$
MW : 262.31

TEST	SPECIFICATION
Appearance	Red to rusty brown powder
Dye content	$\geq 95 \%$
Melting Point	147 – 157° C
Loss on drying	$\leq 0.5 \%$
Residue on Ignition	$\leq 0.5 \%$
Insoluble in Toluene	$\leq 1.0 \%$
Chloride (Cl)	$\leq 0.1 \%$
Solubility: (1 mg/ml Chloroform)	Clear to slightly hazy, Orange to dark red to dark red-brown solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

