



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Flavianic Acid Hydrate
Catalog# : S10086
CAS : 483-84-1
CI number : -
Formula : $C_{10}H_6N_2O_8S \cdot xH_2O$
MW : 314.23 (anhydrous basis)

TEST	SPECIFICATION
Appearance	Yellow powder or crystalline powder
Assay: (titration with NaOH)	96.5% - 103.5% (anhydrous basis)
Melting Point	148° C - 151 ° C
Loss on drying	< 15.0%
Solubility: (0.1% in H ₂ O)	Clear to slightly hazy, yellow to yellow-orange solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

