



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Gallocyanine
Catalog# : S10093
CAS : 1562-85-2
CI number : 51030
Formula : $C_{15}H_{13}ClN_2O_5$
MW : 336.73

TEST	SPECIFICATION
Appearance	Black powder
Dye content	~90.0%
Lambda max.	620nm in ethanol
Solubility: (10 mg/ml, 1N NH4OH)	Opaque dark purple solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

