



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Guinea Green B
Catalog# : S10096
CAS : 4680-78-8
CI number : 42085
Formula : $C_{37}H_{36}N_2O_6S_2$
MW : 690.82

TEST	SPECIFICATION
Appearance	Dark Green to Maroon-Purple Powder
Lambda Max	618 - 620 nm in water
Solubility:	1 % in H ₂ O; Dark Blue to Dark Blue-Green Solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

