



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Indigo Synthetic
Catalog# : S10098
CAS : 482-89-3
CI number : 73000
Formula : $C_{16}H_{10}N_2O_2$
MW : 262.26

TEST	SPECIFICATION
Appearance	Dark Blue to Violet to Purple Powder
Identification (IR)	Conforms to standard
Loss on drying	$\leq 1 \%$
Dye Content	$\geq 90 \%$

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

