



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Janus Black
Catalog# : S10101
CAS : 4443-99-6
CI number : 11825
Formula : $C_{28}H_{26}ClN_5O$
MW : 483.99

TEST	SPECIFICATION
Description	Brown to Black powder
Dye content	> 70%
Lambda max.	584 nm
Melting Point	200° C

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

