



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Lacmoid
Catalog# : S10104
CAS : 42249-61-6
CI number : -
Formula : $C_{24}H_{16}N_2O_6$
MW : 428.39

TEST	SPECIFICATION
Appearance	Dark blue, dark violet to black powder
Lambda max	611nm in methanol
Solubility: (1 mg/ml MeOH)	Opaque, dark blue solution
Visual transition interval	From pH 4.4 (red) To pH 6.4 (blue)

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

