



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Leishman Stain  
Catalog# : S10105  
CAS : 12627-53-1  
CI number : -  
Formula : -  
MW : 774.96

TEST	SPECIFICATION
Appearance	Dark-greenish black powder
Solubility: (0.1% MeOH)	Clear blue solution with greenish fluorescence
Absorptivity 1%/ 1cm	
at Lambda max. 525.0 nm	> 650
at Lambda max. 641.0 nm	> 900
Loss on drying	< 8%
Suitability for Microscopy	Passes Test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

