



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Light Green S.F.Certified
Catalog# : S10107
CAS : 5141-20-8
CI number : 42095
Formula : $C_{37}H_{34}N_2O_9S_3Na_2$
MW : 792.86

TEST	SPECIFICATION
Appearance	Deep Purple Powder
Lambda Max	630 (422) nm (in Water)
Melting Point	288 ° C
Dye Content	≥ 80 %

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

