



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Martius Yellow
Catalog# : S10112
CAS : 605-69-6
CI number : 10315
Formula : $C_{10}H_6N_2O_5$
MW : 234.17

TEST	SPECIFICATION
Appearance	Orange to brownish-yellow to yellow powder
Loss on drying	$\leq 4.5\%$
Solubility: (1 mg/ml methanol)	Clear, yellow to dark yellow solution
UV-Vis (C = 0.01 g/l in Methanol)	E (370 \pm 3) nm = 6,000 minimum
	E (293 \pm 3) nm = 9,000 minimum
	E (249 \pm 3) nm = 14,000 minimum

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

