



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Methylene Violet 3RAX  
Catalog# : S10120  
CAS : 4569-86-2  
CI number : 50206  
Formula :  $C_{22}H_{23}ClN_4$   
MW : 378.91

TEST	SPECIFICATION
Appearance	Green crystalline Powder
Melting Point	285 ° C
Lambda Max	557 nm in Water

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

