



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Nigrosin/Alcohol Soluble
Catalog# : S10130
CAS : 11099-03-9
CI number : 50415
Formula : -
MW : -

TEST	SPECIFICATION
Appearance	Black Powder
Identification (IR)	To conform to structure
Absorption max (In Ethanol)	565 ± 3nm

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

