



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Nigrosin W/S  
Catalog# : S10131  
CAS : 8005-03-6  
CI number : 50420  
Formula : -  
MW : -

TEST	SPECIFICATION
Appearance	Black Crystalline Powder
Lambda max.	570 nm (In 50 % Ethanol)
Solubility: 1 mg/ml, H <sub>2</sub> O	Opaque black with purple cast solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

