



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Nile Blue A Certified
Catalog# : S10133
CAS : 3625-57-8
CI number : 51180
Formula : $C_{40}H_{40}N_6O_6S$
MW : 732.84

TEST	SPECIFICATION
Description	Green to greenish-black powder
Melting Point	> 300 deg. C (dec.)
Dye content	>65%
Lambda max	633 nm (in 50 % Ethanol)
Solubility	
H ₂ O	20 mg/ml
EGME	8 mg/ml
EtOH	3 mg/ml

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

