



TriStains

**DAWN SCIENTIFIC**

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

## Product Specification

Product : Nuclear Fast Red  
Catalog# : S10136  
CAS : 6409-77-4  
CI number : 60760  
Formula :  $C_{14}H_8NO_7SNa$   
MW : 357.28

TEST	SPECIFICATION
Appearance	Dark-red to Red-brown powder
Identification (FTIR)	To conform to standard
Assay (UV vis): 10mg/ltr. Water (N/100 soda Ash solution)	At 535.0+/- 2nm = 0.406
Solubility: (4 mg + 4 ml of H <sub>2</sub> O)	Clear dark red solution
Suitable for Determination of Calcium	To Pass Test

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

