



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Oil Red O Certified
Catalog# : S10137
CAS : 1320-06-5
CI number : 26125
Formula : $C_{26}H_{24}N_4O$
MW : 408.51

TEST	SPECIFICATION
Description	Dark Greenish-Brown to Black to bluish-red Powder
Identification	To pass test
Solubility 1mg/1ml of Ethanol:Chloroform (1:1)	Clear dark red solution
Dye content	≥ 70.0 %
Loss on drying (at 105 Deg. C/1 hr.)	< 5.0 %
Lambda Max	518 (359) nm (in Toluene)
Melting Point	120 ° C

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

