



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Para Red
Catalog# : S10145
CAS : 6410-10-2
CI number : 12070
Formula : $C_{16}H_{11}N_3O_3$
MW : 293.28

TEST	SPECIFICATION
Description	Red powder
Dye Content	~95 %
Loss on drying	< 5%
Melting Point	248 - 252 deg. C
Solubility 0.1% in toluene	Clear to slightly hazy red to orange solution
Lambda max.	488nm in toluene

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

