



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : Patent Blue V
Catalog# : S10148
CAS : -
CI number : 42051
Formula : $C_{27}H_{31}N_2O_7S_2Na$
MW : 582.7

TEST	SPECIFICATION
Description	Dark blue powder
Dye content	> 85.0 %
Subsidiary dyes	< 2.0%
Water-Insoluble Matter	< 0.2%
Leuco Base	< 4.0%
Combined Ether Extract	< 0.2%
Arsenic	< 3 mg/kg
Heavy metals (as Pb)	< 40 mg/kg
Lead (Pb)	< 10 mg/kg
Mercury (Hg)	< 1 mg/kg

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

