



TriStains

DAWN SCIENTIFIC

YOUR SCIENTIFIC PARTNER

ISO 9001:2015 CERTIFIED

Product Specification

Product : 4-Phenylazophenol
Catalog# : S10150
CAS : 1689-82-3
CI number : 11800
Formula : $C_{12}H_{10}N_2O$
MW : 198.2

TEST	SPECIFICATION
Appearance	Yellow to gold powder
Assay	$\geq 98\%$
Lambda max.	347nm in methanol
Melting Point	150 °C to 152 °C (lit.)
Solubility: 1 mg/ml. EtOH	Clear, dark yellow solution

- For laboratory use only.
- The contents of Product specifications are subject to change.
- Storage: For information refer respective product safety data sheet (SDS).

TriStains

